SimulatorTM
Owner's manuel

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PHYSIPRO

This booklet contains important information about this product.
Please give to the owner upon delivery.

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Physipro Inc. is proud to count you among its customers and would like to thank you for the confidence you have shown by purchasing our product.

This owner's manual was created to provide you with all the information needed to allow you to use our product in a safe and optimal way. For all adjustments and settings, we strongly recommend you contact your distributor.

For Physipro Inc. your satisfaction is our highest priority.

Product information
Purchase date:
Serial number:
Distributor:
Address:
Telephone :

The product presented in this manual may not be an exact representation of the product you have on hand. The instructions in this manual are still valid for your product.

Physipro reserves the right to modify or change this manual without notice.

Table of content

Description	4
Components	4
Technical specifications	5
Recommendations	6
Warnings	8
Instructions	9
Wheel locks	9
Tilt angle	10
Backrest	10
Adjusting the backrest angle.....	10
Gridded backrest.....	11
Tension adjustable backrest.....	13
Rigid backrest for Foam-in-Place Seating (FIPS).....	14
Inserting postural components.....	14
Seat	15
Adjusting seat depth.....	15
Adjusting the pelvic supports.....	15
Anti-thrust wedge measurement system.....	16
Abductor mechanism and pommel.....	16
Headrest.....	17
Armrest	19
Footrest	20
Optional tools.....	22
Cleaning and Maintenance.....	23
Cleaning	23
Maintenance	25
Maintenance checklist.....	26
Storage	28
Warranty.....	29
Français.....	31

ENGLISH

Description

This assessment tool is used to determine the exact manufacturing criteria required for a postural support devices in the aim of meeting the specific needs of an individual. Characterized by its remarkable versatility, the Simulator allows the health care professional to try out different adjustments as well determine which postural components, be it for support or correction, are best suited for their client.

Equipped with integrated measurement tools, the Simulator facilitates taking an individual's anthropometric measurements and makes it easier to determine the correct tilt and backrest inclination angles needed to achieve proper seat angle orientation in relation to the ground.

Most of the Simulator's adjustments are performed quickly without the need of any tools, this makes the device an indispensable assessment tool.

Components

TECHNICAL SPECIFICATIONS	
Maximum weight capacity	265 lb (120 kg)
Tilt mechanism	Adjustable dynamic tilt angle: From 0° to 30°
Reclining Gas Spring Backrest	Adjustable backrest inclination angle: From 85° to 165°
Gridded Backrest	Can be removed to allow the installation of a tension adjustable backrest or a rigid backrest for Foam-in Place Seating (FIPS)
Lateral supports	Adjustable height : 11" (28 cm) to 18" (46 cm) Distance between laterals: From 8 ½" (22 cm) to 20" (51 cm)
Rigid seat	Adjustable seat depth: From 14" (36 cm) to 22" (56 cm)
Pelvic supports	Adjustable width: From 7" (18 cm) to 20" (51 cm)
Anti-thrust wedge measuring system	Adjustable thickness: From 0" (0 cm) to 2 ½" (6.4 cm)
Abductor mechanism and pommel	Pivoting and removable Height and depth adjustments
Symbio multi-axe headrest support	Adjustable height: 8" (20 cm) adjustment range Laterally adjustable: 6" (15 cm) adjustment range - 3" (8 cm) on the left and right Adjustable depth: From -4" (-10 cm) to 5" (14 cm)
Armrest	Adjustable height: From 0" (0 cm) to 16" (40 cm) Adjustable width: From 8" (20 cm) to 21" (53 cm) Adjustable depth: 10" (25 cm) adjustment range
Footrest	Adjustable height: From 12" (30 cm) to 21" (53 cm) Adjustable width: From 12" (30 cm) to 30" (76 cm) Adjustable angle: 60°, 70° and 90°
Optional Components and Accessories	Premium tension adjustable backrest Rigid backrest for Foam-in Place Seating (FIPS) Lateral supports (Straight, Contour or Extended) Postural components for seat cushions and backrests Lumbar strap Lap tray
Pediatric Components and Accessories	Abduction pommel Lateral control headrest with integrated & adjustable extensions Two-point padded positioning belt Footplates Premium tension adjustable backrest Lateral supports (Straight, Contour or Extended) Lumbar strap Lap tray

Recommendations

WARNING Do not use this equipment without having read and understood this owner's manual in its entirety. This booklet contains essential information and instructions to ensure the safety of the client and any other individual who is near the Simulator.

Safety Inspection Checks

A thorough inspection of the Simulator should be completed upon delivery, as well as, at the intervals indicated in the Maintenance Checklist, to ensure the client's safety.

- ✓ Verify upon delivery, that the Simulator is assembled with all the components selected in the order form and that no parts are missing or damaged.
- ✓ Verify that the Simulator rolls smoothly and functions normally.
- ✓ Make sure all parts operate without abnormal noises, vibrations, or irregular movements. The presence of one of these conditions may indicate that a part is not adequately fastened or that the Simulator has sustained damage.
- ✓ Make sure the propulsion wheels and caster wheels are firmly attached and that no object interferes with their functional performance.
- ✓ Verify that wheel locks function effectively.
- ✓ Make sure anti-tips are securely attached and function as intended.
- ✓ Verify that both the seat and the backrest are securely attached to the frame.
- ✓ Make sure armrests are securely fastened and locked in place.
- ✓ Make sure footrests are securely fastened.
- ✓ Make sure the belt is firmly attached to the frame.

Safety Guidelines

Several safety measures must be taken to ensure the safety of the occupant and the individuals who are near the Simulator. The following list is not exhaustive. It is the responsibility of every individual to remain cautious in all actions undertaken.

- ✓ The Simulator is an assessment tool and is not designed to be pushed over long distance.
- ✓ When transferring your client from their wheelchair to the Simulator, ensure that wheel locks are engaged.
- ✓ Never attempt to surmount a tall obstacle, changes in height will affect the stability of the Simulator.
- ✓ Never sit more than one person at a time in the Simulator.
- ✓ Do not exceed the weight limit of the Simulator. Exceeding the weight limit increases the risk of losing control of the Simulator, which could result in serious injury to the user or others as well as damage to the Simulator chair. (To verify the maximum load of the Simulator, please refer to the Technical Specifications section of this manual;
- ✓ Never attach objects to the backrest, except those provided with the Simulator.
- ✓ The client must never attempt to tilt the Simulator without assistance.
- ✓ The client must never use the footrest as a platform during a transfer or when attempting to lift oneself.
- ✓ Never use detachable or removable parts to lift the Simulator. Always use the rigid parts of the frame when lifting.
- ✓ Be careful not to pinch your fingers when installing a component or using a removable part.
- ✓ Always be aware of your environment. Taking a few moments to familiarize yourself with your surroundings, before using the Simulator in a new environment, will allow you to avoid obstacles and potential dangers.
- ✓ Avoid exposure to excessive moisture. Always dry the Simulator if it gets wet or when you have finished cleaning it;
- ✓ Avoid exposing the Simulator to extreme temperatures to prevent the risk of hypothermia or burns. The Simulator's surfaces can become frigid in cold temperatures and extremely hot in high temperatures or if exposed to direct sunlight.

ENGLISH

Warnings

Warnings for healthcare professionals

- Physipro recommends that you follow the instructions provided in this manual and familiarize yourself with the Simulator's operating instructions before beginning a client's assessment.
- Physipro recommends that the positioning belt be buckled when a client is seated in the Simulator.
- Always engage wheel locks when an individual is left alone in the Simulator.
- Always maintain a good posture when lifting or tilting the Simulator. Remember to keep your back straight and to bend your knees when lifting.


WARNING

Failure to observe these warnings and precautions increases the risk of losing control of the Simulator, which may result in serious injury to the client or other individuals and may damage the assessment wheelchair.

Transfers

To safely transfer client from one surface to another, the following procedure should be followed.

1. Position the Simulator as close as possible to the transfer surface, with the front caster wheels pointed forward.
2. Immobilize the Simulator by engaging wheel locks.
3. Remove the footrest and ask the client to place their feet on the ground.
4. Remove or swing away the armrest that is located on the side of the intended transfer.
5. Transfer the client.

Note - The use of a transfer board is recommended to ensure the user's safety.

WARNINGS

- Be aware that during the transfer there is a moment when the Simulator is no longer underneath the client.
- The client should never use the footrest as a platform during a transfer or when attempting to lift themselves, as this may cause the Simulator to tip forward, potentially resulting in serious injury and damage to the chair.

Instructions

Required tools

The Simulator is equipped with multiple adjustment handles to facilitate the various modifications that can be made to the assessment wheelchair, however a 4 mm Allen key will be required for some adjustments.

Wheel locks

Using wheel locks

To lock the wheels

Push both handles **A**, located on each side of the wheelchair, forwards.

To unlock the wheels

Pull both handles **A**, located on each side of the wheelchair, backwards.

Physipro recommends that wheel locks be engaged during an assessment.


ENGLISH


Tilt angle

The following steps explain how to tilt the assessment wheelchair:

1. Press and hold the activation levers **A** located under the push handles.
2. Increase or decrease tilt angle.
3. Release the levers **A** at the desired tilt angle.

NOTE: You can read the tilt angle on the goniometer, located on the lower part of the backpost.

Backrest

Adjusting the backrest angle

The following steps explain how to adjust the backrest angle:

1. Press and hold both levers **B** located under the tension bar.
2. Increase or decrease the backrest angle.
3. Release levers **B** when the desired angle is obtained.


NOTE: You can read the backrest angle on the goniometer, located under the push handles.

Gridded backrest

The gridded backrest's positioning guide allows you to quickly and efficiently measure the width of an individual's shoulders and chest, as well as determine the correct position of lateral supports.

Adjusting the lateral supports

Lateral supports are positioned by default in the outer adjustment slots. Before beginning adjustments, determine if the lateral supports should be placed in the inner adjustment slots. If so, follow these steps :

1. Unscrew both handles **A** completely.
2. Move both lateral supports to the inner adjustments slots.
3. Refasten handles **A**.

Once the lateral supports are positioned in the correct adjustment slots, you may begin adjusting the lateral supports:


1. Loosen both handles **A**.
2. Adjust the height of the lateral supports by moving the laterals upwards or downwards along the adjustment slots.
3. Adjust the distance between both lateral supports.
4. Firmly tighten handles **A** when lateral supports are correctly positioned.
5. Note the position of the laterals on the gridded backrest.

ENGLISH

Replacing the backrest

It is possible to remove the gridded backrest to install a tension adjustable backrest or a rigid backrest for Foam-in-Place Seating (FIPS) during the assessment. To remove the gridded backrest follow these steps:


1. Press triggers **A** and secure them in hooks **B**, this action will release the spring loaded latch that holds the upper backrest anchors **E** in the slots of the backrest supports **C**.
2. Pull the backrest up so that the lower backrest anchors **F** are completely removed from hooks **D**.

To install an adjustable tension backrest or a rigid backrest for Foam-in-Place Seating (FIPS), follow these steps:

Note: To install the selected backrest, the triggers **A** on the backrest supports **C** must be secured into hooks **B**.

1. Insert the lower backrest anchors **F** into hooks **D**.
2. Raise the backrest slightly to insert the upper backrest anchors **E** into the slots of the backrest supports **C**.
3. Press the backrest downwards so that both spring latches snap in place.

Tension adjustable backrest

The tension adjustable backrest is available in adult or pediatric sizes. This backrest can be adjusted to meet the needs of the client and also allows the installation of lateral supports and postural components.

Adjusting backrest tension


To adjust the tension of the backrest you must tighten or loosen each strap individually according to the morphology of the client's back.

WARNING Do not tighten straps excessively, this may deform the backrest.

Adjusting the lateral supports

Height adjustment

1. Loosen handle **A**.
2. Position the lateral support at the desired height.
3. Tighten handle **A**.

Depth adjustment

1. Loosen screws **B** with a 4 mm Allen key.
2. Position the lateral support at the desired depth.
3. Firmly tighten screws **B**.

Width adjustment

1. Loosen handle **C**.
2. Turn the mechanism until the desired width is obtained;
3. Tighten handle **C**.
4. Loosen handle **A** to adjust the angle of the lateral support and then retighten the handle.

ENGLISH**Rigid backrest for Foam-in-Place Seating (FIPS)**

The installation of a rigid backrest for Foam-in-Place Seating (FIPS) on the Simulator facilitates the use of a molding bag system. The openings in the backrest allows the healthcare professional to pass their hands through the backrest to ensure a successful impression of their client's back is made.

To install the rigid backrest for Foam-in-Place Seating (FIPS) please refer to the steps outlined in the section: *Replacing the backrest* on page 12.

To use the Foam-in-place seating system, please refer to the instruction manuel.

**Inserting postural components**

Postural components, such as wedges and sacrolumbar correctors, can be added to the backrest during the assessment.


Seat

Adjusting seat depth

To adjust seat depth follow these steps:

1. Press and hold handle **A**.
2. Pull or push the seat until the desired depth is obtained.
3. Release handle **A**.

Note: A ruler is located on the side of the seat to facilitate seat depth measurement.


Adjusting the pelvic supports

Pelvic supports are adjustable in width and in depth.

Width adjustment

1. Loosen handles **B**.
2. Position the pelvic supports as desired;
3. Tighten handles **B**.

Depth adjustment

1. Loosen handles **C**;
2. Adjust the depth of the pelvic supports;
3. Tighten handles **C**.


ENGLISH

Anti-thrust wedge measurement system

This system helps determine the adequate thickness of an anti-thrust wedge to limit forward migration.

1. Turn handle **D**, located at the back of the Simulator, until the desired height is obtained;
2. With a ruler, measure the distance between the seat plate and the top of the footrest support block.


Abductor mechanism and pommel

The push button abductor mechanism can be adjusted in height and moved laterally for proper positioning of the pommel. This mechanism can also pivot under the seat when the front button is pressed.

Adjusting the height of the pommel

1. Loosen handle **E**.
2. Adjust the height of the pommel.
3. Tighten handle **E**.

Adjusting the lateral position of the pommel

1. Loosen handle **F** located under the seat.
2. Slide the mechanism to the left or right until properly positioned;
3. Tighten handle **F**.

Pivoting the pommel under the seat

1. Press button **G**.
2. Pivot the pommel under the seat.


Headrest

The Symbio headrest support can be moved laterally on the tension bar and offers height, depth and angle adjustments to properly position the headrest. The extensions of lateral control headrest can be adjusted in depth and at an angle.

Adjusting the support's lateral position

1. Loosen all four (4) screws **A** with a 4 mm Allen key.
2. Slide the headrest support to the left or right until properly positioned.
3. Tighten all four (4) screws **A**.


Adjusting the height of the headrest

1. Loosen handle **B**.
2. Adjust the height of the headrest.
3. Tighten handle **B**.

Adjusting the depth of the headrest

1. Loosen handles **C** and **D**.
2. Adjust the depth of the headrest so it supports the clients head.
3. Tighten handles **C** and **D**.

Note: It is likely that headrest height will need to be readjusted.

Adjusting the headrest's lateral position


1. Loosen handle **E**;
2. Move the headrest to the left or right until properly positioned.
3. Tighten handle **E**.

ENGLISH



Adjusting the angle of the headrest

1. Loosen all three (3) screws **F** with a 4 mm Allen key.
2. Move the headrest until properly angled.
3. Tighten all three (3) screws **F**.

Adjusting the lateral extensions

1. Loosen both screws **G** and both screws **H**.
2. Adjust the position and depth of the lateral extension.
3. Tighten all screws **G** and **H**.
4. Loosen screw **I**.
5. Adjust the angle of the lateral extension.
6. Tighten screw **I**.
7. Repeat steps for the other side.

Headband

The headband is offered as an option. It is installed with velcro fasteners that are attached to the lateral extensions.


Note: It is important that the headband is properly adjusted. The headband should be tight enough to prevent the client's head from dropping forward but should not cause discomfort. Headband should be used under supervision of a caregiver.

Armrest

The Simulator's armrest are adjustable in height, width, depth and angle.


Adjusting the height and width

1. Unlock lever **A**.
2. Raise or lower the armrest until the desired height is reached.
3. Turn the armrest inwards or outwards to adjust the width between both armrests.
4. Lock lever **A**.

Adjusting the depth and angle

1. Loosen handle **B**;
2. Move the armrest until the desired depth and angle is obtained.
3. Tighten handle **B**.

Lap tray

The lap tray is offered as an option. It is installed with velcro fasteners that are attached to the armrests, as shown in the photo.


ENGLISH

Footrest

The footrests are adjustable in height, width and angle, they can also be removed to facilitate transfers. The footplates are adjustable in depth and angle.


Adjusting the height of the footrest

1. Press and hold lever **A**.
2. Move the footplate upwards or downwards until the desired height is obtained.
3. Release lever **A**.

Adjusting the width

1. Loosen handle **B** on both footrests.
2. Move footrests to the left or right until the desired width between both footrest is obtained.
3. Tighten both handles **B**.


Removing the footrests

1. To remove the footrests, unscrew and take off end cap **C** from tube.
2. Loosen handle **B**.
3. Pull the footrest outwards until completely removed from the tube.

Adjusting the footrest angle

1. Pull on the release pin **D**;
2. Adjust the angle of the footrest by reinserting the release pin **D** in the desired adjustment hole.

Adjusting the footplates

1. Loosen handle **E**.
2. Position the footplate at the desired depth and adjust its angle to properly support the clients foot.
3. Tighten handle **E**.

ENGLISH

Optional assessment tools

Physipro offers additional measuring tools and evaluation kits that can be used with the Simulator during a postural assessment.

Vernier caliper

Accurate and precise, the vernier caliper is ideal for taking anthropometric measurements.

Available in two sizes:
25" (63 cm) and 50" (130 cm)

Foam kit

This foam kit allows the healthcare professional to try out different combinations of foams in the aim of identifying the best cushion composition required to provide maximum support and comfort to their client.

Magnetic goniometer

The magnetic goniometer can be used to measure different postural angles.

Evaluation bag

The evaluation bag comes with various postural components of different dimensions for seat and backrest, that allow the healthcare professional to evaluate which components are best suited for their client.

Cleaning and Maintenance

Regular cleaning and maintenance will extend the lifespan of the Simulator. We strongly recommend that when cleaning your assessment wheelchair, you inspect all the parts and components to ensure they are in good condition.

General Recommendations

- Immediately remove all contaminants (food and/or organic).
- After an assessment, disinfect the seat cushion, armrests, backrest, and any other surface with a disinfectant spray.

Cleaning

Frame and painted surfaces

A weekly cleaning is recommended.

1. Clean painted surfaces with a mild soap or a neutral based detergent and water (2 ounces (6 cl) for 8L of water) and rinse thoroughly.
2. Dry with a clean cloth to absorb excess moisture.

WARNINGS

- Never use an abrasive cleaner.
- Never use a pressure washer.

Axle and Moving Parts

A weekly cleaning is recommended.

1. Remove all dirt and dust from the axles and moving parts.
2. Clean the axles and moving parts with a damp cloth and dry with a clean cloth.

Rigid backrest

A weekly cleaning is recommended.

1. Clean the rigid backrest with a cloth that is slightly dampened with warm water and a mild detergent.
2. Wipe with a clean dry cloth.

ENGLISH**Tension adjustable backrest, cushion, lateral supports and headrest**

A weekly cleaning is recommended.

Cleaning the outer covers

1. Remove cover and fasten the zipper.
2. Machine wash in cold water on the gentle cycle or hand wash with a mild detergent.
3. Allow to air dry completely.
4. Reinstall cover.

**WARNING**

Do not iron.

Cleaning foam padding

1. Remove cover.
2. Wipe foam padding lightly with a damp cloth, without soap or detergent.
3. Use a clean dry cloth to gently press on padding to absorb excess moisture, then let air dry completely.
4. Reinstall cover.

**WARNINGS**

- Do not use soap.
- Do not submerge in water.
- Never machine wash or dry the foam padding.

Comfort Accessories

A weekly cleaning is recommended.

1. Gently wipe with a cloth dampened with a mild soap or neutral detergent.
2. Remove excess moisture by lightly patting the surfaces with a clean cloth.
3. Allow to air dry completely.

Disinfection

Disinfection helps destroy infectious agents and eliminate pathogenic micro-organisms with the use of bactericidal or virucidal cleaning products.

1. Wipe surfaces with disinfecting wipes containing at least 70% alcohol.
2. Let product react for 15 minutes.
3. Rinse with a clean, water-dampened cloth.
4. Remove excess moisture by lightly patting the surfaces with a clean cloth.
5. Allow to air dry completely.

WARNING

Never use a bleaching agent on fabrics or straps.

Maintenance

Performing regular maintenance will extend the Simulator's lifespan and ensure your safety, as well as your clients, during use. The Simulator's Maintenance checklist can be found on the following pages.

WARNINGS

- Upon receiving the Simulator, a thorough inspection of all components and parts should be completed to ensure their proper function.
- A regular inspection of the Simulator's components and parts, at the intervals indicated in the Maintenance Checklist, is necessary to keep your assessment wheelchair in good condition.
- The inspection by a qualified technician is recommended twice a year to ensure the Simulator's proper function.

ENGLISH

Maintenance Checklist	At reception	Each week	Each month	Every 6 months
General				
■ Check the solidity of the Simulator's frame and components	✓			✓
■ Make sure all parts operate without abnormal noises, vibrations, or irregular movements.	✓	✓		
■ Make sure the Simulator rolls in a straight line	✓	✓		
■ Check tires for wear or deformation	✓		✓	
■ Make sure wheels roll smoothly and caster forks rotate properly	✓		✓	
■ Check cushion, covers and comfort upholstery for wear	✓			✓
■ Check the effectiveness and condition of wheel locks	✓		✓	
Tilt mechanism				
■ Make sure activation levers and tilt mechanism functions properly	✓		✓	
■ Make sure gas cylinders are not leaking oil	✓		✓	
■ Check condition of cables and make sure they are well attached	✓		✓	
Backrest				
■ Make sure levers and the backrest inclination mechanism functions properly	✓		✓	
■ Make sure the lateral supports adjustment handles function properly	✓		✓	
■ Make sure the backrest is securely attached	✓		✓	
■ Make sure the gridded backrest is not damaged	✓			✓

Maintenance Checklist	At reception	Each week	Each month	Every 6 months
Seat				
■ Make sure the seat depth mechanism functions properly	✓	✓		
■ Make sure the pelvic support adjustment mechanism functions properly	✓	✓		
■ Make sure the anti-thrust wedge measurement system functions properly	✓	✓		
■ Make sure the abductor mechanism functions properly	✓	✓		
Headrest				
■ Check the effectiveness of the headrest support handles	✓	✓		
■ Check the effectiveness of the headrest support adjustment mechanisms	✓	✓		
Armrest				
■ Check the effectiveness of the armrest supports handles and levers	✓	✓		
■ Check the effectiveness of the armrest supports adjustment mechanisms	✓	✓		
Footrest				
■ Check the effectiveness of the footrest supports handles and levers	✓	✓		
■ Check the effectiveness of the footrest supports adjustment mechanisms	✓	✓		
Two-point padded positioning belt				
■ Check the condition and make sure the positioning belt is securely attached	✓		✓	
■ Make sure the belt buckle functions properly	✓	✓		

ENGLISH**Repair Services**

If replacement parts are needed, please contact our customer service at 1 800 823-2252 or at order@physipro.com for the complete list of authorized repair facilities and the repair procedure.

Certain parts and components, such as the wheels, armrests, footrests, headrest and seat cushion can be removed and sent directly to a repair facility, after receiving a written approval by Physipro Inc.

Certain distributors may provide replacement units during the repair process. For more information, contact your representative.

Physipro Inc. will provide repair parts for a minimum of 5 years or will provide a compatible replacement option.

Replacement Parts

For more information, you can contact us by phone at 1 800-668-2252, by email at order@physipro.com or visit us in person. Physipro Inc. is located at 370, 10th Avenue South, Sherbrooke (Quebec) Canada J1G 2R7.

Storage

The Simulator must be stored in a clean and dry area.

**WARNINGS**

- Never store the Simulator in excessively humid or hot conditions or in direct sunlight.
- Never store the Simulator near sharp objects or in an area where it may be struck by an other object.
- If the Simulator has been stored for several weeks, please consult the Maintenance Checklist, and proceed to a thorough inspection before use. It is important to verify the proper functioning of each component.
- If the Simulator has been stored for over two (2) months, an inspection by a qualified technician is required.

Disposal of Medical Devices

Do not burn or throw away medical devices outside of designated areas.

Life span

The estimated life span of the Simulator is five (5) years.

WARRANTY

Physipro Inc. is committed to maintain all devices, components, and parts in good working condition, and to repair all design or manufacturing defects, that may impede on the proper functioning or breakdown of the product, at its expense.

This warranty applies to all devices, components and parts and covers all materials and workmanship required for the complete execution of repair. This warranty also includes administrative and transport costs incurred for repair services.

The warranty periods for the parts and components of the Simulator and the optional assessment tools are as stated:

- 5 years** - Frame and cross-brace against defects in materials and workmanship.
- 12 months** - Parts and labor for original defects in material and workmanship of the Simulator's components and optional components.
- 3 months** - Vernier caliper, magnetic goniometer, foam kit and evaluation bag
- 30 days** - Tires, comfort upholsteries, armrests, seat and backrest covers, and ball bearings.

Devices and components that are replaced or repaired under one of these warranties will remain subject to the warranty for the remaining warranty period.

Replacement and repairs carried out during the warranty period must be made with the original parts and components.

Physipro inc. agrees to repair or replace any inoperative or defective part for the duration of the warranty period.

The warranty does not cover damages:

- Attributed to the installation of a component or part, provided by a third party or a manufacturer, without the supplier's authorization.
- Caused by a third party or a manufacturer, that does not have the supplier's authorization.
- Arising from improper use or failure to comply to the operation and maintenance instructions and the recommendations provided in the owner's manual.
- Damages that occur during transportation.

ENGLISH**Warranty service**

If a defect is detected during the warranty period, please contact the retailer from whom you purchased your Physipro Inc. product as soon possible. Please refrain from using your product until an inspection has been made by a qualified technician.

In the event you do not receive satisfactory warranty service, please write directly to Physipro Inc. at the address listed below or send us an e-mail at **order@physipro.com**, please include the retailer's name and address, a brief description of the defect, the invoice date, and the products serial number, if applicable. DO NOT RETURN product prior to receiving a written return authorization from Physipro Inc.

The defective unit or parts must be returned for warranty inspection within thirty (30) days of the return authorization date at the following address:

Physipro Inc.
370, 10th Avenue South, Sherbrooke, Quebec, J1G 2R7

Exclusions and limitations

Compliance with the safety precautions and recommendations mentioned in this manual is crucial to avoid injury to the user or other persons. Inadequate or inappropriate adjustments will increase the risk of accident. It is important to thoroughly read the manual before making any adjustments. Physipro Inc. will not be held responsible for bodily injury or property damage caused by misuse or incorrect adjustments made to its product by an unqualified person. Physipro Inc. disclaims all liability for damages that occur during transportation or are the result of situations beyond its control.

The foregoing warranty is exclusive and in lieu of all other expressed warranties.

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