## **CHEST SUPPORT WITH CENTER BUCKLE**

### **Owner's Manual**

This booklet contains important information about this product.  
Please give to the owner upon delivery.

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**ENGLISH****INTRODUCTION**

Physipro Inc. is proud to count you among its customers and would like to thank you for the confidence you have shown by purchasing our product.

This owner's manual was created to provide you with all the information needed to allow you to use our product in a safe and optimal way. For all adjustments and settings, we strongly recommend you contact your distributor.

Physipro Inc., relieves itself of all liability should physical injury or property damage occur due to the lack of maintenance or misuse of our product, or from modifications made to the product without the prior written consent of Physipro Inc.

**PRODUCT INFORMATION**

**Date of purchase:** \_\_\_\_\_

**Supplier:** \_\_\_\_\_

**Address:** \_\_\_\_\_

**Telephone:** \_\_\_\_\_

The product presented in this manual may not be an exact representation of your product. However, the instructions provided in this manual are still valid.

Physipro Inc. reserves the right to modify or change this manual without notice.

**TABLE OF CONTENT**

|                                   |           |
|-----------------------------------|-----------|
| DESCRIPTION.....                  | 3         |
| COMPOSITION.....                  | 3         |
| TECHNICAL SPECIFICATIONS.....     | 4         |
| RECOMMENDATIONS.....              | 5         |
| WARNINGS .....                    | 5         |
| INSTALLATION .....                | 6         |
| Installing the lower straps.....  | 6         |
| Installing the upper straps.....  | 7         |
| Adjusting the chest support ..... | 8         |
| MAINTENANCE.....                  | 9         |
| WARRANTY.....                     | 10        |
| <b>FRANÇAIS .....</b>             | <b>11</b> |

## DESCRIPTION

The chest support with center buckle provides stability and facilitates the process of removing the occupant from their wheelchair. Offered in three different styles, each adapted for either men, women, or children, these chest supports are ergonomically designed to comfortably hold the occupant's torso upright when seated, without hindering their upper body movement. Made with a soft neoprene padding with a lycra trim, the chest support features multiple slide buckles for accurate adjustment as well as, four easy-to-release plastic buckles to facilitate transfers.

This item meets the requirements of California 116 & 117, NFPA 701, CE marking and NF EN1021-1, EN1021-2 norms.

## COMPOSITION


Neoprene  
Lycra  
Nylon straps  
Auto buckle  
Plastic buckles  
Plastic slide buckles  
Strap anchor plates  
Cam buckles

## Installation hardware options

Belt attachment blocks for round seat rails  
Rigid extension plates

## ENGLISH

## TECHNICAL SPECIFICATIONS

| Size          | Chest Support - Men |        |        |        |
|---------------|---------------------|--------|--------|--------|
|               | A                   | B      | C      | D      |
| <b>Small</b>  | 234 mm              | 376 mm | 181 mm | 123 mm |
| SC7730        | 9 ¼"                | 14 ¾"  | 7"     | 4 ¾"   |
| <b>Medium</b> | 276 mm              | 467 mm | 203 mm | 112 mm |
| SC7731        | 10 ¾"               | 18 ¼"  | 8"     | 4 ½"   |
| <b>Large</b>  | 330 mm              | 502 mm | 203 mm | 95 mm  |
| SC7732        | 13"                 | 19 ¾"  | 8"     | 3 ¾"   |

| Size          | Chest Support - Women |        |        |        |
|---------------|-----------------------|--------|--------|--------|
|               | A                     | B      | C      | D      |
| <b>Small</b>  | 321 mm                | 381 mm | 173 mm | 133 mm |
| SC7736        | 8"                    | 15"    | 6 ¾"   | 5 ¼"   |
| <b>Medium</b> | 343 mm                | 422 mm | 203 mm | 133 mm |
| SC7737        | 13 ½"                 | 16 ½"  | 8"     | 5 ¼"   |
| <b>Large</b>  | 483 mm                | 476 mm | 440 mm | 121 mm |
| SC7738        | 19"                   | 18 ¾"  | 6 ½"   | 4 ¾"   |

| Size            | Chest Support - Children |        |        |        |
|-----------------|--------------------------|--------|--------|--------|
|                 | A                        | B      | C      | D      |
| <b>One Size</b> | 170 mm                   | 280 mm | 215 mm | 130 mm |
| SC7739          | 6 ¾"                     | 11"    | 8 ½"   | 5"     |

## RECOMMENDATIONS

**To ensure the longevity of your chest support, the following precautions should be applied.**

- Before use, it is important to verify the quality and condition of the chest support's components.
- Physipro Inc. recommends that a health care professional such as a doctor or therapist, experienced in seating and positioning, be consulted to install, and adjust the chest support.
- Clean regularly and inspect the chest support's components periodically.
- Remove any foreign object (crumbs, debris, etc.) from the chest support.

**It is recommended to consult your doctor or therapist in the month following the installation of a new postural support device to make sure it is efficient.**

## WARNINGS

- Physipro Inc. chest supports are reserved exclusively for use on mobility aids (manual or power wheelchairs).
- Always use the provided mounting hardware to install your chest support.
- Never push or lift wheelchair with the chest support, doing so may cause the chest support to detach from the wheelchair which may result in damage or injury. Always lift wheelchair with the rigid parts of the frame.
- The chest support should be tested before prolonged use. To avoid redness due to pressure points, it is necessary to regularly check for redness in the various critical pressure areas, especially around bony prominences.
- Before use, let the chest support get to room temperature to limit discomfort caused by excessively hot or cold surfaces.
- Avoid exposing the chest support to extreme cold or heat and avoid prolonged contact to fluids or moisture. Never place near open flame.
- Crash tests have not been conducted on this product and even if the mobility aid in which it is installed has obtained the ISO-717619, it will not be maintained.

## ENGLISH

### INSTALLATION


The chest support can be installed on several types of wheelchairs. Depending on your wheelchair's seat rail type (round or square tubing), follow the proper installation procedure for the lower anchor plates, the upper anchor plates can be installed directly to a rigid backrest or to a tension bar. Follow the installation procedure of your choice.

#### Installing the lower straps

##### Round seat rails

**Required tools:** 4mm and 5 mm Allen keys

1. Determine the correct location of the anchor plates on the seat rails.
2. Position both parts of the belt attachment block to the seat rail.
3. With a 5mm Allen key, tighten screw **B** firmly.
4. With a 4mm Allen key, secure the anchor plate to the attachment block by fastening screw **A**.
5. Repeat these steps for the other side.


**Note** – Make sure both anchor plates are placed in the corresponding adjustment holes of the opposite seat rail.

##### Square seat rails

**Required tool:** 4mm Allen key

1. Determine the correct location of the anchor plates on the seat rails.
2. Tighten screw **A** with a 4mm Allen key.
3. Repeat these steps for the other side.


**Note** – Make sure both anchor plates are placed in the corresponding adjustment holes of the opposite seat rail.

## Installing the upper straps

### Rigid backrest (without extension plate)

**Required tool:** 3mm Allen key

1. Drill four ¼" holes in the rigid backrest.
2. Insert sleeves **A** into the holes.
3. Position the cam buckles **B** directly on the rigid backrest, align the cam buckles with the holes.
4. Insert screws **C** into sleeves **A**.
5. Tighten screws **C** firmly.
6. Insert upper straps into the cam buckles **B**.


### Rigid Backrest (with extension plate)

If your wheelchair has a low backrest, an extension plate can be installed to reduce pressure on your shoulders.

**Required tool:** 3mm Allen key

1. Drill four ¼" holes in the rigid backrest.
2. Insert sleeves **A** into the holes.
3. Align the extension plates with the holes.
4. Insert screws **C** into sleeves **A**.
5. Tighten screws **C** firmly.
6. Align the cam buckles **B** with the upper holes of the extension plates.
7. Insert screws **D** through the holes of the extension plates, through the cam buckles **B**.
8. Tighten screws **D** firmly with nuts **E**.
9. Insert upper straps into the cam buckles **B**.


## ENGLISH

### Tension bar


**Required tools:** 4mm and 5mm Allen keys

1. Position both parts of the attachment block on the tension bar
2. Insert screws **A** and **B** into the attachment block.
3. With a 5mm Allen key, tighten screw **B** firmly.
4. Position the adjustment plate **C** at the desired height and tighten screws **A** with a 4mm Allen key.
5. Repeat these steps for the other attachment block.

### Adjusting the chest support

To adjust the chest support, you must slide the four adjustable straps through the plastic buckles to lengthen or shorten the straps.

#### WARNINGS

- When using your wheelchair, it is important to always fasten your chest support for added protection.
- It is important to adjust your chest support to maximize comfort and safety.
- Improper use can result in serious injury or death.

## MAINTENANCE

Regular cleaning eliminates the accumulation of dirt and grime and helps extend the life of your chest support. A weekly cleaning is recommended.

### Cleaning the neoprene padding

1. Remove the neoprene padding from the straps by unfastening all four plastic buckles.
2. With a mix of water and a soft soap or neutral detergent, gently wipe the neoprene padding and then rinse by wiping with water only.
3. Remove excess moisture by lightly patting the neoprene padding with a clean cloth and then let air dry completely.

#### WARNINGS

- Do not submerge in water, this may damage to auto-buckle
- Never machine wash or dry.
- Do not iron.

### Cleaning the straps

1. With a mix of water and a soft soap or neutral detergent, gently wipe the straps.
2. Remove excess moisture by lightly patting the straps with a clean cloth.
3. Let air dry completely.

#### WARNINGS

- Never machine wash or dry.
- Do not iron.

### Cleaning the hardware

Regularly dust off hardware with a clean cloth.

### Storage

Postural supports must be stored in a dry and secure area.

#### WARNINGS

- Never store your postural support in direct sunlight, near a heat source or in a high-humidity environment.
- Never place your postural support near sharp objects or in an area where it may be struck by an object.

### Life span

The estimated life span of your postural support is three (3) years.

### Disposal of Medical Devices

Do not burn or throw away medical devices outside of designated areas.

**ENGLISH****WARRANTY**

The warranty period for the parts and components of the chest support with center buckle is **12 months**.

This limited warranty covers defects in material and workmanship on all the standard or optional parts and components of this product, from the date of shipment, for the duration stated above.

Parts and components that are replaced or repaired under one of these warranties will remain subject to the warranty for the remainder of the warranty period. Physipro Inc. agrees to repair or replace, free of charge, any defective part or component for the full duration of the warranty period.

**This warranty does not apply to:**

- a) Products that were misused, neglected or have not been maintained and stored in accordance with the recommendations of the manufacturer.
- b) Products that were intentionally or accidentally altered, or products that have been damaged through repairs or modifications made without the prior written consent of Physipro Inc. or of an authorized retailer.
- c) Normal wear and tear or damages caused by the user such as punctures, tears, or burns.

**Warranty service**

If a defect is detected during the warranty period, please contact the retailer from whom you purchased your Physipro Inc. product as soon possible. Please refrain from using your product until an inspection has been made by a qualified technician.

In the event you do not receive satisfactory warranty service, please write directly to Physipro Inc. at the address listed below or send us an e-mail at **order@physipro.com**, please include the retailer's name and address, a brief description of the defect, the invoice date, and the products serial number, if applicable. **DO NOT RETURN** product prior to receiving a written return authorization from Physipro Inc.

The defective unit or parts must be returned for warranty inspection within thirty (30) days of the return authorization date at the following address:

**Physipro Inc.**

370, 10th Avenue South, Sherbrooke, Quebec, J1G 2R7

**Exclusions and limitations**

Compliance with the safety precautions and recommendations mentioned in this manual is crucial to avoid injury to the user or other persons. Inadequate or inappropriate adjustments will increase the risk of accident. It is important to thoroughly read the manual before making any adjustments. Physipro Inc. will not be held responsible for bodily injury or property damage caused by misuse or incorrect adjustments made to its product by an unqualified person.

Physipro Inc. disclaims all liability for damages that occur during transportation or are the result of situations beyond its control.

**The foregoing warranty is exclusive and in lieu of all other expressed warranties.**

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